





# 260-BURKIT02 Canine Bur Set, FG



All burs are 19 mm and used in High-speed handpieces. Suited for treatments on dogs.




Round burs -To be used for caries excavation, cavity preparation, root canal access and osseous surgery.

Burs	#	Quantity	Description
	260-H131014	2	Round, size 4, 19 mm
	260-H131018	2	Round, size 6, 19 mm


Cross-cut taper fissure burs – Ideal for sectioning multi-rooted teeth and reducing crown height.


Burs	#	Quantity	Description
	260-H33L31012	1	Flat-End Cross-Cut Fissure, # 701L, 19 mm
	260-H33L31016	2	Flat-End Cross-Cut Fissure, # 702L, 19 mm

Round-End Fissure – For sectioning teeth and removing broken roots in gingiva.

Burs	#	Quantity	Description
	260-H21R31012	1	Round-End Fissure, # 1158, 19 mm

Inverted cone – With a slightly rounded edge on blade corners, which produces a smoother and superior cut. Generally used for entering carious pits, producing an undercut in cavity preparations with a flat floor.

Burs	#	Quantity	Description
	260-H231010	1	Inverted cone, # 35, 19 mm

Burs	#	Quantity	Description
	260-886231012	1	Flame Shape short, diamond (fine grit), size 5, 19 mm



D I S P O M E D