






260-BURKIT03 Feline Bur Set, FG



The lengths in this set vary from 16 to 19 mm to cover the general use in the treatment of cats and small dogs.




Round burs -To be used for caries excavation, cavity preparation, root canal access and osseous surgery.

Burs	#	Quantity	Description
	260-H133005	1	Round, size 1/4, 16 mm
	260-H133014	2	Round, size 4, 16 mm
	260-H133023	1	Round, size 8, 16 mm


Cross-cut taper fissure burs – Ideal for sectioning multi-rooted teeth and reducing crown height.


Burs	#	Quantity	Description
	260-H3333012	2	Flat End Taper Cross-Cut Fissure, # 701, 16 mm
	260-H33L31012	1	Flat-End Cross-Cut Fissure, # 701L, 19 mm

Round-End Fissure – For sectioning teeth and removing broken roots in gingiva.

Burs	#	Quantity	Description
	260-H33R31016	1	Round End Taper Cross-Cut Fissure, # 1702, 19 mm

Inverted cone – With a slightly rounded edge on blade corners, which produces a smoother and superior cut. Generally used for entering carious pits, producing an undercut in cavity preparations with a flat floor.

Burs	#	Quantity	Description
	260-H233006	1	Inverted cone, # 33 1/2, 16 mm

Burs	#	Quantity	Description
	260-36831016	1	Flame Shape short, diamond (fine grit), size 5, 19 mm



D I S P O M E D