CERTAINTY WHERE IT MATTERS MOST






Capnostream[™] 35 Portable Respiratory Monitor

CONTINUOUS MONITORING SIMPLE AND EASY

The Capnostream[™] 35 respiratory monitor is designed with features to assist clinicians to help enhance efficient workflow.

- Portable, rugged, lightweight, ergonomic design helps you facilitate moving a patient between areas of care
- Smart algorithms engineered to give you an early indication of changes in respiratory status that may need intervention
- Continuous monitoring of etCO_2, respiration rate, and \mbox{SpO}_2





COMPREHENSIVE MONITORING ACROSS THE CONTINUUM OF CARE

Quick, easy, and reliable access to respiratory status. It's the key to help provide safe and effective care from a patient under sedation during a procedure to monitoring a patient on the general care floor. With a combination of Microstream[™] capnography and Nellcor[™] pulse oximetry technologies, you can count on the Capnostream[™] 35 respiratory monitor to deliver. Continuous monitoring of etCO₂, respiration rate, and SpO₂. In virtually any clinical setting. It brings certainty where it matters most.

o 11				
General characteristics				
Unit dimensions	213 mm (h) x 137 mm (w) x 55 mm (d) (8.38 in (h) x 5.39 in (w) x 2.16 in (d))			
Unit weight	1 kg. 2.2 lb.			
Power supply				
Input power	100-240 VAC, 50/60 Hz, ± 10%			
Battery				
Battery type	Lithium ion			
Battery pack nominal voltage and energy	7.2V; 18.72Wh			
Battery operation	Three hours for removable battery pack; 20 minutes for internal battery pack			
Battery charging time	Up to five hours when the monitor is off and up to eight hours when the monitor is on, for both batteries			
Microstream [™] capnography				
CO ₂ , etCO ₂ range	0-150 mmHg			
CO2 accuracy	0-38 mmHg: ± 2 mmHg 39-150 mmHg: ± (5% of reading + 0.08 for every 1 mmHg above 38 mmHg)			
Respiration rate range	0-150 bpm			
Flow rate	50 (42.5 \leq flow \leq 65) ml/minute, flow measured by volume			
Nellcor [™] pulse oximetry				
SpO₂ measurement range	1-100%			
Pulse rate range	20 to 250 beats per minute			

PM35-VET-KIT Contents:

Optional Accesories:

Item Code	Description	Qty.	Item Code	Description
PM35MN02	Capnostream [™] 35 Portable Respiratory Monitor	1	CS03743	Pole Clamp
V-SAT	Vet Sensor	1	PM35VAA	Cap 35 vesa adaptor assembly
DOC10	Pulse Oximetry Cable, DOC10	1	D-YS	Dura-y sensor (For sensor replacement)

VISIT CAPNOSTREAM35.COM* FOR MORE INFORMATION

medtronic.com*

*This is a global website. It is not specific to Canada.

Claims made throughout this brochure are based on clinical trials. © 2016 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. Currently licensed under Covidien LLC. Manufactured by: Covidien LLC, 15 Hampshire Street, Mansfield, MA, 02048 USA. CA-AH-0010-E Rev. 2017/01

Medtronic