

4.2. OVERALL DIMENSIONS

The following figures give the overall dimensions of the possible configurations allowed, depending by the configuration that has been ordered:

4.2.1.FRONT VIEW (REST POSITION) - BOTTOM MOUNT




Fig. 4.2

A	
40cm (16'') bracket	79cm (31'')
80cm (31'') bracket	119cm (46'')
110cm (43'') bracket	149cm (59'')

4.2.2.FRONT VIEW (REST POSITION) - TOP MOUNT




Fig.4.3

B	
40cm (16") bracket	79cm (31")
80cm (31") bracket	119cm (46")
110cm (43") bracket	149cm (59")

4.2.3.SIDE VIEW (OPEN) - TOP MOUNT




Fig.4.4

C1 - LONG CONE	
40cm (16") bracket	140cm (55")
80cm (31") bracket	180cm (71")
110cm (43") bracket	210cm (83")
C2 - SHORT CONE	
40cm (16") bracket	150cm (59")
80cm (31") bracket	190cm (75")
110cm (43") bracket	220cm (87")
C3 - FULL SIZE	
40cm (16") bracket	188cm (74")
80cm (31") bracket	228cm (90")
110cm (43") bracket	258cm (102")

4.2.4.SIDE VIEW (CLOSED) - TOP MOUNT




Fig.4.5

D	
40cm (16'') bracket	60cm (24'')
80cm (31'') bracket	100cm (39'')
110cm (43'') bracket	130 (51'')

4.2.5.SIDE VIEW (OPEN) - BOTTOM MOUNT




Fig.4.6

E1 - LONG CONE	
40cm (16") bracket	140cm (55")
80cm (31") bracket	180cm (71")
110cm (43") bracket	210cm (83")
E2 - SHORT CONE	
40cm (16") bracket	150cm (59")
80cm (31") bracket	190cm (75")
110cm (43") bracket	220cm (87")
E3 - FULL SIZE	
40cm (16") bracket	188cm (74")
80cm (31") bracket	228cm (90")
110cm (43") bracket	258cm (102")

4.2.6.SIDE VIEW (CLOSED) - BOTTOM MOUNT




Fig.4.7

F	
40cm (16'') bracket	60cm (24'')
80cm (31'') bracket	100cm (39'')
110cm (43'') bracket	130 (51'')